LEAD LAB:

Community Based Exploratory Learning in the Context of Teaching Instrumentation With a Single Analyte

Alanah Fitch
Loyola University Chicago
afitch@luc.edu
8/01

Lead Lab: Teaching Instrumentation with a Single Analyte

"[The suspected material] is to be acted on with dilute nitric acid, and the resulting liquid filtered. If the addition of dilute sulphuric acid to the filtered liquid produce a white, if chromate of potash a yellow, and, lastly, if sulphuretted hydrogen a black precipitate, the presence of lead is fully determined."

Mitchell, 1848 A.D.: An Early Analytical Text

Take two to three pounds of the juice of the Lunaria, add a solution of lead, alum, salpetre, sal ammoniac, silver-litharge, quick-silver sublimate, vinegar and ginger, then distil and calcinate. This powder is then thickened with a silver solution into paste. One part in a hundred of this paste, heated on a copper plate, could theoretically turn tin, lead, or quicksilver into silver at a ratio of one to a hundred parts.

Raymond Lully, 1235 A.D.

Contents

A.	Table of	of Contents
	Key Co	oncept Sequences
	Class S	Syllabusv
		Attached Article Teamwork, C&EN, Nov. 15, 1999
B.	Selecti	ng A Method of Lead Analysis
A. B. C.		ments
٠.	1:	Population Analysis
		Attached Article H. L. Needleman et. al. The New Eng. Journal of Medicine, 1979, 30, 13, 689.
	2:	S/N and filtering
	3:	Lead by lead sulphide (turbidity)
	4:	Lead by lead chromate (gravimetric)
	5:	Commercial spot tests
	6*:	Electrolytic deposition of PbO ₂ (gravimetric)
	7:	Anodic Stripping Voltammetry (alloy formation, electrochemical detection) (EPA Method 7472, 7063) 47
	7.	Attached Article: T. Liu, D. Lai, J. D. Osterioh. Anal. Chem., 1997, 69, 3539-3543
		Attached Article M. Zelic and M. Lovric, Electrochem. Acta, 35, 11/12, 1990, 1701-1706
		Attached Article B. J. Feldman, J. D. Osterlog, B. H. Hata, and A. D'Alessandro, Anal. Chem., 1994, 66,
	0	1983-1987.
	8:	PbS ion selective electrode (potentiometric, based on metal sulphide alloy)
	0.4	Attached Article: S.Kamata and K. Onoyama, Anal. Chem., 1991, 63, 1295-1298
	9*:	Titrametric, EDTA with Eriochrome Indicator
	10*:	Titrametric, EDTA, UV-Vis detection
	11:	Lead Iodide determination
	12:	Atomic Absorption, Graphite Furnace (EPA 200.9)
		Attached Article J. Creed, T. Martin, L. Lobring, and J. O'Dell, Env. Sci. Tech., 1992, 26, 102-106.
		Attached Article B. J. Marquardt, S. R. Goode, J. M. Angel, Anal. Chem., 1996, 68, 977-981
	13:	Chelation by Dithizone, UV-Vis determination (ASTM E40-77)
		Attached Material: Excel for Chemists: Deconvolution
		Attached Article H.M.N.H. Irving, D. C. Rupainwar, and S. S. Sahota, Anal. Chim. Acta, 1969, 45, 249-
		254.
	14:	Stability of Dithizone by Cyclic Voltammetry
		Attached Article L. Tomcsanyi, Anal. Chim. Acta, 1977, 88, 371-376.
		Attached Article A. Y. Kassim, Inorg. Chim. Acta, 1978, 27, 243.
	15.	Fluorescence Quenching of Calcein Blue by Lead
	16	Fluorescence of Eyrothrocyte protoporphyrin
		Attached Article: P. B. Hammond, R. L. Bornschein, and P. Succop. Env. Ressearch, 1985, 38, 187-196.
	17:	IR study of Lead-EDTA chelate
	18:	NMR of ²⁰⁷ Pb-EDTA complex
	19:	ICP-MS determination of isotopes of lead (EPA 200.8)
	20*:	Preparative Chromatography for Lead84
		Attached Article: W. I. Manton, Chem. Geology, 1988, 73, 147-152.
	21:	Ion Pair HPLC of Lead EDTA
		Attached Article; Y. Shijo, K. Takada, N. Uehara, Anal. Sci., 1993, 6, 315.
		Attached Article W. Liu, H. K. Lee, J. of Chromatography A, 1998, 796, 385-395
	22:	GC of tetraethyl Lead87
		Attached Article M. R. Brunelto, J. L. Burguera, M. Burguera, and D. Chakraborti, Atom. Spectroscop.,
		1992, 13, 4, 123.
	23:	Hot Plate Acid Digestions (EPA 200.2)89
	24*:	Sequential Soil Extractions 91

	25:	Soil Sampling	93						
	26:	House Dust Sampling: Wet Wipe	95						
	27*:	House Dust Sampling: Vacuum							
	28:	Water Sampling	98						
	29*:	Blood Sampling	99						
	30*.	Hand Dust Sampling	100						
D.	Readii	ngs in Ethics	101						
	31.	Community Based Problem Solving							
	32.	The Needleman Papers	105						
	33.	The Johnson Controls Case	106						
		Attached Material: Adapted from Markkula Center for Applied Ethics, Santa Clara University							
		Attached Artices:							
		C. B. Ernhart, S. Scarr, D. F. Geneson, Ethics and Behavior, 1993, 3, 1, 73-93							
		H. M. Needleman, Ethics and Behavior, 1993, 3, 1, 95-101							
		S. Scarr and C. B. Ernhart, Ethics and Behavior, 1993, 3, 2, 199-206							
		H. L. Needleman et. al., New England J. of Medicine, 1990, 322,2, 83							
		U.S. District Court. F. D. Wisonsin #84-C-472							

E. Appendices and Useful Information

U.S. Supreme Court 499 US 187

^{*} not yet student proofed

Key Concept Sequences

The following experiments have been tried and tested over a 6 year lab sequence for undergraduate chemistry majors. They are arranged in a historical and conceptual sequence allowing students to acquire a background understanding of chemistry.

Concept Sequence 1:

Sampling statistics control the **overall** analytical procedure beginning with in field sampling both in randomizing individual samples and in collecting a sufficient data base of information to derive meaningful information. Sampling statistics also control the **limit of detection**, a figure of merit used to select an instrument for use in analysis. By increasing the number of instrumental measurements using computer based collection of information, the standard deviation of the measurement goes down, and, therefore, so does the limit of detection.

Students are asked to make a running set of comparisons of the instrumental methods throughout the semester, so they need to be able to determine the **limit of detection** on each instrument and understand what factors come into play. **Experiments 1 and 2** set the student up to properly understand statistics and to use computer spread sheets, as well as some simple digital and electronic filtering to reduce noise.

Concept Sequence 2:

The final analysis is dependent upon the weakest link. Often that weak link is the analyst's understanding of the chemistry that precedes the instrumentation. One recurring theme in lead analysis is the need to control the pH in order to avoid lead hydroxide formation. Students are asked to build on the spread sheet skills in producing an alpha plot for lead hydroxides and to compare that data to measured free lead vs pH (ISE, experiment 8). Chemistry continues to be introduced in the form of chelation and phase separation in the UV-Vis Dithizone experiment. If the pH is not properly controlled the experiment can not work. The chemistry in this experiment also depends upon selectivity in chelation and on phase separation, concepts that reappear in the IR, NMR, and Chromatographic experiments.

Concept Sequence 3:

Lead exists in a variety of isotopes. The NMR and MS labs exploit the fingerprint pattern of lead.

Concept Sequence 4:

Ethics is implicitly and explicitly built into this lab. **Explicit** ethical issues to be explored include "professional" ethical issues relating to professional relationships within the work groups, honesty, integrity (Lab 1: a good lab book; Lab 23: digestion and designing internal quality control). We also explore how statistics (Lab 1) informs public policy in a different manner than it does analytical measurements. We directly explore (last lab) how a scientific disagreement should be carried out. **Implicit** ethic issues concern ownership of data between unequal collaborators (community groups and scientists), the extent to which scientists have a moral obligation to impart scientific knowledge to the general public, and the obligation that scientists have to not participate in scare science and/or to "use" community groups as mere learning tools.



Part A: Class Syllabus and Information

Instrumental Methods of Analysis Chemistry 310/311 Fall 2001 Loyola University Chicago.

Goal: To be able to chose an instrument for lead analysis using limits of detection, and, reliability, and ease of use. Problem Based Learning to monitor house dust with collaborating 6th grade children or to produce a publication quality research project on lead around a municipal waste incinerator.

Things students will be graded on:

Exams (4), all must be taken, however lowest grade may be averaged with highest to get a value for the fourth exam . 400 points

Informed Class participation (20%) 100 points

A: 90%

B 80- 89.9% C. 70-79.9% D. 60-69.9%

Text Teaching Instrumentation with a Single Analyte (Fitch)

Optional: Skoog, Leary, Nieman 6th Ed.

EXAM FORMAT

3 essays, 4 problem solving Problems derived from homework and from lab Essay may include reading a contemporary article

List of Equations on the following pages may be xeroxed, highlighted, **but not annotated**, and brought to exams.

CHEAT SHEETS

This may be xeroxed & brought to exams. You may use a highlighter on this sheet only.

Constants

 $h = 6.62608 \times 10^{-34} \text{ J-s}$

 $k = 1.38 \times 10^{-23} \text{ J/K}$

 $c = 3x10^8 \text{ m/s}$

R = 8.31451 J/mol-K

 $F = 9.6485 \times 10^4 \text{ C/mol}$

 $N = 6.022 \times 10^{23}$ particles/mole

The Golden Key: The Normal Equation

$$f(x) = \frac{1}{s\sqrt{2p}}e^{-\frac{1}{2}\left(\frac{x-m}{s}\right)^2}$$

Electronics and Noise and Statistics

$$V = q/C$$

dq/dt = I

 $V_o = -\hat{a}(\mathring{a}_{\scriptscriptstyle{-}} + \mathring{a}_{\scriptscriptstyle{+}})$

$$X_c = \frac{1}{2pfC}$$

$$Z = \sqrt{R^2 + X_c^2}$$

$$\Delta f = \frac{1}{t_{r}}$$

$$R = \frac{\left(\overline{x}_a - \overline{x}_b\right)}{6s}$$

$$\frac{S}{N} = \frac{\overline{x}_a - \overline{x}_b}{s}$$

$$LOD = x_{blank} \pm 3s_{blank}$$

$$LOQ = x_{blank} \pm 9s_{blank}$$

$$6s \approx \left| Peak_{max} \right| - \left| Peak_{min} \right|$$

Optics

$$nl = c$$

$$E = h\mathbf{n}$$

$$h_i = \frac{c}{n}$$

$$f = \frac{\mathbf{n}_m 2\mathbf{n}}{c}$$

$$I = \frac{2t\mathbf{h}}{n}$$

$$f = F/d$$

$$W = \frac{\Delta \mathbf{l}}{2D^{-1}} \qquad \mathbf{l}_{\text{max}} = \frac{2898}{T}$$

$$n\mathbf{I} = d(\sin i + \sin r)$$

$$J_{I} = \left(\frac{2\mathbf{p}hc^{2}}{\mathbf{I}^{5}}\right) \left(\frac{d_{I}}{e^{\frac{hc}{kTI}} - 1}\right)$$

$$P\mathbf{a}\left[\frac{1}{f/\#}\right]^2 \qquad R = nN$$

$$R = nN$$

$$\frac{\sin \boldsymbol{q}_1}{\sin \boldsymbol{q}_2} = \frac{\boldsymbol{h}_2}{\boldsymbol{h}_1}$$

$$\frac{I_R}{I_o} = \left(\frac{\boldsymbol{h}_2 - \boldsymbol{h}_1}{\boldsymbol{h}_2 + \boldsymbol{h}_1}\right)^2$$

$$I_s = \left[\frac{I_o}{\boldsymbol{l}^4}\right] \left[\frac{8\boldsymbol{p}^4\boldsymbol{a}^2}{r^2}\right] \left[1 + \cos^2 \boldsymbol{q}\right]$$

Spectroscopy (100-800 nm)

Spectroscopy
$$> 1\mu$$

$$A = -\log T = -\log\left(\frac{P}{P_o}\right) = \log P_o - \log P = \mathbf{e}bC$$

$$A = \log \left[\frac{P_{o, \mathbf{I}_1} + P_{o, \mathbf{I}_2}}{P_{o, \mathbf{I}_1} (10^{-\mathbf{e}_1 bC}) + P_{o, \mathbf{I}_2} (10^{-\mathbf{e}_2 bC})} \right] \qquad \mathbf{m} = \frac{m_1 m_2}{m_1 + m_2}$$

$$m = \frac{m_1 m_2}{m_1 + m_2}$$

$$\frac{s_c}{C} = \frac{0.434 s_T}{T \log T} \qquad \left| \Delta \mathbf{I}_{1/2} \right| = \frac{\mathbf{I}^2 \Delta \mathbf{n}}{c}$$

$$\mathbf{n} = \frac{1}{1} = \frac{1}{2pc} \sqrt{\frac{k}{m}} = 5.3x10^{-12} \sqrt{\frac{k}{m}}$$

$$\Delta \mathbf{I} = 2 \left(\frac{\mathbf{I}}{c}\right) \left(\frac{2kT}{\mathbf{p}m}\right)^{1/2} = 4.848x10^{-7} \mathbf{I} \sqrt{\frac{T}{m}}$$

#bands possible: 3N-5 or 3N-6

$$\frac{N_j}{N_c} = \left(\frac{P_j}{P_c}\right) e^{\frac{-\Delta E}{kT}}$$

$$\Delta N = 2b(\boldsymbol{n}_1 - \boldsymbol{n}_2)$$

$$b = \frac{\Delta N}{2(\boldsymbol{n}_1 - \boldsymbol{n}_2)}$$

Electrochemistry

$$\Delta G = -nFE = -RT \ln K$$

$$E = E^{o} - \frac{0.0591}{n} \log \frac{[\text{Re } d]}{[Ox]}$$

$$E_{ISE} = const' + 0.0592 \log(a_{analyte})$$

$$E_{measured} = E_{cell} = E_{rhs} - E_{lhs} = E_{cat} - E_{an}$$
 $i_c = AC_d \left[\frac{\P dE}{\P t} \right] = AC_d \mathbf{n}$

$$\frac{hc}{I} \propto \Delta E^{o}$$

$$i = nFAD \frac{\P C}{\P x} = nFAD \frac{\Delta C}{d}$$

where
$$d \approx \sqrt{2Dt}$$

$$I_{\text{limiting},RDE} = 0.620 nFAD^{2/3} \mathbf{n}^{-1/6} w^{1/2} C$$

$$I_{polarography} = 0.607 nm^{2/3} D^{1/2} C_{bulk} t^{1/6}$$

$$i_c = AC_d \left| \frac{\P dE}{\P t} \right| = AC_d \mathbf{n}$$

$$q = nFN$$

The Diffusion Equation: Does it Look Familiar?

$$C_{x,t} = \frac{C_o}{\sqrt{4pDt}} e^{-\left(\frac{x^2}{4Dt}\right)}$$
$$= \frac{C_o}{\sqrt{2Dt}\sqrt{2p}} e^{-\frac{1}{2}\left(\frac{x}{\sqrt{2Dt}}\right)^2}$$

Chromatography

$$k' = \frac{t_r - t_m}{t_m}$$

$$k' = \frac{KV_s}{V_M}$$

$$t_{r,B} - t_m$$

$$K$$

$$\mathbf{a} = \frac{t_{r,B} - t_m}{t_{r,A} - t_m} \qquad \mathbf{a} = \frac{K_B}{K_A}$$

$$N = \frac{16t_r^2}{W^2} \qquad H = \frac{L}{N}$$

$$H = \frac{W^2 L}{16t_r^2}$$

$$H = \frac{B}{u} + C_s u + C_m u$$

$$B = 2yD_m$$

$$C_s = \frac{qk'd_f^2}{(1+k')^2 D_s}$$

$$C_m = \frac{f(d_p^2, d_c^2, u)}{D_m}$$

$$R_s = \left[\sqrt{\frac{N}{4}}\right] \left[\frac{\mathbf{a} - 1}{\mathbf{a}}\right] \left[\frac{k'_b}{k'_b + 1}\right]$$

$$t_{r,b} = \left[\frac{H}{u}\right] 16R_s^2 \left[\frac{\mathbf{a}}{\mathbf{a} - 1}\right]^2 \left[\frac{(1+k'_b)^3}{k'_b^2}\right]$$

Schedule: Alanah Fitch, Fall, 2001									
Day	M	T	W	Th	F				
8:30	Chem 311	Chem 311	Chem 311	Babysitting					
9:30	Lab	Lab	Lab		At Env.				
10:30					Studies/Sci. Program				
11:30				Chem. Seminar	Damen 111 phone: 88992				
12:30	Lunch	Lunch	Lunch	Lunch	Office Hours				
1:30	At Env.		Chemistry						
2:30	Studies/Sci. Program	Lecture	Research	Lecture	Chem 101				
3:30	Damen 111	Chem 310		Chem 310					
4:30	office Hours								

Office Department Phone email

402 Flanner Hall Chemistry 508-3119 afitch@luc.edu

111 Damen Hall Environmental Studies/Sci. 508-8992

Home 847-256-8243 (between 6-9 p.m.)

Office Hours: M 1:30 to 5 F 9:30-2

Revised Schedules: Fall 2000 Instrumental Analysis/ September 21, 2000									
						Stat	Che m	Isoto pe	Ethics
1	Aug 29	Т	Society Experiment 1 Divide into Semester						
	Aug 31	Th	Statistics and Ethics	Teams I and II Alanah	Teams I and II Alanah				
2	Sept 5	Т	Signals and Noise	Electronics: Analog and Digital Filtering Experiment 2 Alanah and Carla					
	Sept 7	Th	Signals & Noise, LOD						
3	Sept 12	Т	Relevant Chemistry	Lead Hydroxide Equilibrium Chemistry and Ion Selective Electrodes Experiment 8 Jorge					
	Sept 14	Th	Exam I						
4	Sept 21	Th	Intro Spectros	IR	UV-Vis				
5	Sept 26	Т	Molecular	UV-Vis	IR				
	Sept 28	Th	Spectroscopy	Experiment 13 Experiment 17 Carla Audrius					
6	Oct 3	Т	Fluorescence	Fluorescence	FAA				
	Oct 5	Th	AA	Experiment 15 Carla	Experiment 12 Jorge				
7	Oct 10	Т	AA/Vibrational	FAA Fluorescence Experiment 15 Carla					
	Oct 12	Th	Vibrational Take Home Exam II						
8	Oct 17	Т	Mid Semester Break	Mid semester Break	Mid semester Break				
	Oct 19	Th	NMR						
9	Oct 24	Т	NMR	NMR	ASV				
	Oct 26	Th	Orozco: Demo Sampling	Experiment 18 Audrius	Experiment 7 Jorge				

Oct 27	F		MidTerr	m Grades Are Due							
Oct 31		T		MS	ASV	NMR	10				
Nov 2		Th		Orozco: LUC Field Trip	Experiment 7 Jorge	Audrius	t 18				
Nov 6		M		Last Day to Withdrav	V						
Nov 7		Т		Intro to Electro	Sample Prep						
Nov 9		Th		Orozco: Public Health	Rotation Experiment 23 Rotation Rotation		t 23				
Nov 14		T		Potentiometry	Sample Analysis: Student Decision: Based on Linear Range, LOD, Time, Reliability						
Nov 16		Th		Orozco: Return Data Assist Graphing			7				
Nov 21		Т		ASV/ Take Home Exam III	ICP-MS Experiment 19		t 22				
Nov 23		Th		Thanksgiving	Luke	Audrius					
Nov 28		T		Separations	GC	ICP-MS	. 10				
Nov 30		Th		Separations	Experiment 22 Audrius	Experimen Luke	t 19				
Dec 5		Т		Separations	Ethics & How do we Know th	e "Truth"					
	Oct 31 Nov 2 Nov 6 Nov 7 Nov 9 Nov 14 Nov 16 Nov 21 Nov 23 Nov 28 Nov 30	Oct 31 Nov 2 Nov 6 Nov 7 Nov 9 Nov 14 Nov 16 Nov 21 Nov 23 Nov 28 Nov 30	Oct 31 T Nov 2 Th Nov 6 M Nov 7 T Nov 9 Th Nov 14 T Nov 16 Th Nov 21 T Nov 23 Th Nov 28 T Nov 30 Th	Oct 31 T Nov 2 Th Nov 6 M Nov 7 T Nov 9 Th Nov 14 T Nov 16 Th Nov 16 Th Nov 21 T Nov 23 Th Nov 28 T Nov 30 Th	Oct 31 Nov 2 Th Orozco: LUC Field Trip Nov 6 M Last Day to Withdrav Nov 7 T Intro to Electro Th Orozco: Public Health Nov 14 T Potentiometry Nov 16 Th Orozco: Return Data Assist Graphing Nov 21 T ASV/ Take Home Exam III Nov 23 Th Thanksgiving Nov 28 T Separations Nov 30 Th Separations	Oct 31 Nov 2 Th Orozco: LUC Field Trip Nov 6 M Last Day to Withdraw Nov 7 T Nov 9 Th Orozco: Public Health Nov 14 T Potentiometry Nov 16 Th Orozco: Return Data Assist Graphing Nov 21 T ASV/ Take Home Exam III Nov 23 Th Thanksgiving Nov 28 T Separations GC Experiment 22 Audrius Ethics &	Oct 31 Nov 2 Th Orozco: LUC Field Trip Nov 6 M Last Day to Withdraw Nov 7 T Intro to Electro Nov 9 Th Orozco: Public Health Nov 14 T Potentiometry Nov 16 Th Orozco: Return Data Assist Graphing Nov 21 T ASV/Take Home Exam III Nov 23 Th Thanksgiving Nov 28 T Separations Sample Prep Experiment 23 Rotation Sample Analysis: Student Decision: Interaction and Interaction a	Oct 31	Oct 31	Oct 31	Oct 31 Nov 2 Th Orozco: LUC Field Trip Nov 6 M Last Day to Withdraw Nov 7 T Intro to Electro Nov 9 Th Orozco: Public Health Nov 14 T Nov 16 Th Orozco: Rotation Nov 16 Th Orozco: Rotation Nov 16 Th Orozco: Rotation ICP-MS Experiment 19 Luke Nov 28 T Separations Sample Prep Experiment 23 Rotation Rotation Rotation Orozco: Return Data Assist Graphing CC Experiment 19 Luke ICP-MS Experiment 19 Luke

311 Lab

Text: Lead Lab

Materials required: lab book with carbon pages, 3.5" diskettes

Group work:

Students will work in groups of nor more than 3 and no more than two groups per lab section. The group will be responsible for assigning work associated with the making of standards and the operation of each piece of equipment.

Group composition: Groups will be formed in the first week of lab. Students will interview each other for similarity of commutes and work habits in order to ensure that all members of a group can meet at some mutually agreeable time. Your best friend does not necessarily make your best lab partner in this context.

Scope of Work: Students **must <u>pre-read</u>** the lab. Sometimes work out of assigned class period is required. This will be compensated by cancelling lecture periods.

Each Student Will be monitored to ensure proper rotation of work.

The Weekly Reports

Labor division:

- a) Data collection, tabulation, summarizing
- b) Research
- c) Editing and Copying

are three tasks associated with each report. Each person in the group **must** rotate tasks. Attach the following sheet to each lab (see next).

What the report looks like:

- a) introduction
- b) materials and methods
- C. reduced data (graphs and tables). Graphs are to be inserted immediately after the first reference to the graph in the text. Graphs are to be labeled 1, 2, 3, etc. with a descriptive title.
- d) an essay that is not simply a list of answers to questions
- e) raw data

What you turn in:

- a) The typed report
- b) Work Distribution Form, with division of labor signed by each student

- c) Carbon copy from page of labbook from each student
- d) Individual data as relevant

When you turn it in:

Next lab

Provisional grade given

Edited and redone material with the original submission containing faculty commentary must be re-turned in 1 week after faculty edits.

Two rewrites only are allowed.

What if my group drives me nuts?

The groups can/will be configured three weeks into the semester. Be frank with your colleagues about their effort. If the group doesn't think you are pulling your weight you will end up doing **ALL** the work yourself!!

Lab	Student:	Student:	Student:
Potatos	Data/Graphs		
	Writing/Typing		
	Research/Editing		
ADC	Data/Graphs		
	Writing/Typing		
	Research/Editing		
$Pb(OH)_2$	Data/Graphs		
	Writing/Typing		
	Research/Editing		
Field Sampling	Data/Graphs		
	Writing/Typing		
	Research/Editing		
Sample Prepara			
	Data/Graphs		
	Writing/Typing		
	Research/Editing		
GFAA	Data/Graphs		
	Writing/Typing		
	Research/Editing		
Molecular Abso			
	Data/Graphs		
	Writing/Typing		
	Research/Editing		
IR	Data/Graphs		
	Writing/Typing		
	Research/Editing		
NMR	Data/Graphs		
	Writing/Typing		
M	Research/Editing		
Mass Spec	Data/Graphs		
	Writing/Typing		
AGTI	Research/Editing		
ASV	Data/Graphs		
	Writing/Typing		
CV	Research/Editing		
CV	Data/Graphs		
	Writing/Typing		
IIDI C	Research/Editing		
HPLC	Data/Graphs		
	Writing/Typing		
CC	Research/Editing		
GC	Data/Graphs		
	Writing/Typing Page 27th /Editing		
	Research/Editing		

This page should be inserted with each lab report:

Method	Limit of Detection ug Pb/g solution	Linear Range ugPb/g	Interfences	Ease of Use: Solvent and/or Toxic Materials Used
Gravimetric (Chromate)				
Ion Selective Electrode				
Dithizone				
Calcein Blue				
ASV				
FAA				
GFAA				
ICP-MS				
IR				
NMR				
GC				